

KWIK-STIK™ Plus

Standard Features:

- Two passages from the reference culture
- Canister of 5 KWIK-STIK™ units
- Authentic and traceable
- Shelf-life of 10-24 months
- Refrigerated storage
- Peel-off Certificate of Assay
- Online Certificate of Analysis
- Product warranty



Two Passages From Reference Cultures

Two Passage Microorganisms

MicroBioLogics® offers products designed to enhance your laboratory's quality control program by adding efficiency and reliability. KWIK-STIK™ Plus microorganisms are designed to meet the needs of laboratories requiring QC microorganisms that are a limited number of passages from the reference culture. KWIK-STIK™ Plus microorganisms are two passages from the reference culture.

MicroBioLogics®

217 Osseo Avenue North St. Cloud, MN 56303 Phone: 320.253.1640
800.599.BUGS (2847) Fax: 320.253.6250

www.microbiologics.com



KWIK-STIK™ Plus Product Details

Product Details

KWIK-STIK™ Plus microorganisms have the same characteristics as the KWIK-STIK™ with two additional features:

1. KWIK-STIK™ Plus microorganisms are lyophilized microorganism preparations that are two passages from the reference culture.
2. A peel-off label contains a Certificate of Assay that specifies the identity and traceability of the microorganism to a reference culture, and the number of passages the microorganism preparation has been removed from the reference culture.

KWIK-STIK™ Plus comes in canisters containing five units of a single strain.

Visit www.microbiologics.com or call 320.253.1640 to place your order!